

Patient information – contact lens service

Care of your rigid gas permeable contact lenses

Please read this leaflet carefully as it will instruct you on how best to take care of your rigid gas permeable contact lenses. For further guidance, watch our short film which is available on the Moorfields website:

www.moorfields.nhs.uk/content/contact-lens-insertion-removal-and-care

If there is anything you do not understand in this leaflet, please contact us using the details below:

Phone: 020 7566 2100, Monday-Friday, 9am– 5pm
Email: moorfields.contactlenses@nhs.net
Visit: Clinic 4, contact lens reception desk (Monday to Friday, 9am– 4.45pm).

How can I best take care of my contact lenses?

- Wash and dry your hands thoroughly before inserting or removing your lenses.
- Take care not to catch the lens or your eye with your fingernails – it

is good practice to keep them short.

- Work over a clean, flat surface.
- Use saline solution to rinse your lenses and case.
- To minimise the possibility of mixing your lenses up, it is a good idea to get into the habit of always inserting and removing the same lens first.
- Use the care system or the correct cleaning solutions recommended to you in the clinic.

Please remember not to do any of the following:

- Do not rinse your lenses and case with tap water as this can increase the risk of eye infections.
- Never shower or bath with your lenses in.
- Do not sleep with your lenses in (unless you have been advised to by your contact lens practitioner).
- Avoid swimming with your contact lenses in. This is due to the risk of infection from the water. Prescription swimming goggles and other options may be

available. We advise you discuss this with your contact lens practitioner.

Before inserting your contact lenses

1. Remove the lens from the storage case and make sure it is not damaged, particularly checking the edges for any signs of damage; check too that there are no deposits on the lens (see figure 1 below).
2. Add one drop of the conditioning solution (Boston or Total Care). Gently rub solution over the lens so that both sides are covered.



Figure 1

Inserting your contact lenses

1. Place the contact lens onto the tip of your index finger of your preferred hand.
2. Pull your lower lid down with the middle finger of the same hand.
3. Pull your upper lid up using the middle finger of the other hand. Maintain a firm hold of the upper eyelid to prevent blinking (see figure 2 opposite).

4. Gently place the lens on your eye while firmly holding both eyelids in position. Try to keep both eyes open without blinking.
5. Let go of your **lower** eyelid first and then gently close your **top** eyelid.
6. Once the lens is on your eye, close your eyes slowly and blink a few times.
7. If you wear a lens in the other eye, repeat the steps above.
8. Throw away the solution from your case, rinse it with saline solution and leave it upside down on a clean surface to air dry.

If you have been advised in clinic to use a different method to insert and remove your lenses, ensure you follow those recommended instructions.



Figure 2

Removing your contact lenses

After you remove the contact lens from your eye, it must be cleaned before it is put into the storage case to remove the microbes and deposits that have built up on the lens during the day.

Removal method one

1. Look straight ahead and press the edge of your lower eyelid against the eyeball beneath the contact lens.
2. Lift your top eyelid above the contact lens and then push it down over the eyeball until the eye is shut – the lens should be pushed out and resting between the eyelashes (see figure 3 below).



Figure 3

Removal method two

1. Look straight ahead and open your eye as wide as possible.
2. Place your index finger at the outer corner of your eye and pull your eyelids tight over towards your ear (see figure 4 below).
3. Blink – the lens should come out.



Figure 4

Removal method three

1. If methods 1 and 2 are unsuccessful, a suction holder may be used (see figure 5 on next page).
2. Wet the end of the suction holder with saline and place onto the contact lens. It is imperative that the suction holder is not placed directly onto the surface of the eye.
3. Pull gently to remove the lens from the eye.
4. Remove the suction holder from the lens by sliding it off.



Figure 5

Cleaning your contact lenses

If you have been advised in clinic to use a different method to clean your lenses, ensure you follow those recommended instructions.

1. Place the lens into the palm of your hand and put two or three drops of the **cleaning** solution (Boston, Total Care or Oté cleaner) onto the lens.
2. Rub the lens gently but firmly with the little finger of your

