

Patient information-pharmacy

Compliance aids step by step guides-Autodrop®

The Autodrop® is for use with eye drop bottles. It is designed to help you with the stable positioning of the bottle over your eye, so you can accurately squeeze a drop into your eye.

Instructions for use:

1. After washing your hands, unclip the Autodrop® and place the neck of the bottle into the holder of the clip, with the nozzle facing towards the main body of the aid.
2. Lie down or sit down comfortably, then remove the bottle cap and close the Autodrop® clip. It should clip shut easily.
3. Tilt your head back and lower your bottom eyelid using your finger.
4. Place the aid over your eye with the main body of the aid against your cheek. Look through the pin hole at the top, and when ready, squeeze the bottle to release a drop.
5. Unless your clinician has advised otherwise, close your eye and gently press the inside corner for 30 seconds to help prevent the drop entering your mouth.
6. When finished, open the Autodrop® clip and replace the cap of your eye drop bottle.



Please note: the Autosqueeze® can also be used together with the Autodrop®.

Author: Sarah Thomas, lead satellite pharmacist and Fiona Chiu, associate chief pharmacist
Revision number: 1
Approval date: April 2018

This information can be made available in alternative formats, such as easy read or large print on request. Please call PALS: **020 7566 2325/2324**