

## Patient information-pharmacy

# Compliance aids step by step guides- eye drop bottle opener

The eye drop bottle opener is designed to help open eye bottles, and is ideal for patients who struggle to twist or unscrew eye bottle caps. There are two ends of the opener made in different sizes, which allows the opener to fit and open different sized bottles.

### Instructions for use:

1. Place and gently fit the right sized end of the opener over the bottle cap.
2. Hold the opener in your hand and twist in an anti-clockwise direction to release the cap.
3. To put the cap onto the bottle, hold the opener over the nozzle and twist in a clockwise direction.



Author: Sarah Thomas, lead satellite pharmacist and Fiona Chiu, associate chief pharmacist  
Revision number: 1  
Approval date: April 2018

This information can be made available in alternative formats, such as easy read or large print on request. Please call PALS: **020 7566 2325/2324**